

Anti-CGRP antibody (STJ16101265)

STJ16101265

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CGRP is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-F, IHC-P, IF
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P IF (1:1 000-1:2 000)
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6283
Gene Symbol	S100A12
Uniprot ID	S10AC_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	