

Anti-Adrenomedullin antibody (STJ16101288)

STJ16101288

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Guinea Pig polyclonal antibody anti-Adrenomedullin is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC-F |
| Host/Source | Guinea Pig |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | IHC-F (1:400-1:800) |
| Formulation | Lyophilized; reconstitute in 100 µl dist. water |
| Isotype | |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |