

## Anti-Csk-binding protein antibody (STJ16101301)

STJ16101301

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Csk-binding protein is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	IP, WB, IHC-P
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse, Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Isotype</b>	
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">55824</a>
<b>Gene Symbol</b>	<a href="#">PAG1</a>
<b>Uniprot ID</b>	<a href="#">PHAG1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	