

Anti-C1q antibody {FITC} (STJ16101305)

STJ16101305

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-C1q is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence