

Anti-caspase-1 antibody (C-Term) (STJ16101311)

STJ16101311

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-caspase-1 (C-Term) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC-P
Host/Source	Rabbit
Reactivity	

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	834
Gene Symbol	CASP1
Uniprot ID	CASP1_HUMAN
Immunogen	C-Term
Immunogen Region	C-Term
Specificity	
Immunogen	
Sequence	