

## Anti-caspase-4 antibody (N-Term) (STJ16101312)

STJ16101312

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-caspase-4 (N-Term) is suitable for use in Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	WB, IHC-P
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">837</a>
<b>Gene Symbol</b>	<a href="#">CASP4</a>
<b>Uniprot ID</b>	<a href="#">CASP4_HUMAN</a>
<b>Immunogen</b>	N-Term
<b>Immunogen Region</b>	N-Term
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	