

## Anti-Pan Cytokeratin antibody (STJ16101360)

STJ16101360

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Pan Cytokeratin is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	IHC-P, IHC-F
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Sheep, Pig, Dog, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-P (1:100-1:400) IHC-F
<b>Formulation</b>	TRIS with 0.03% sodium azide, pH7, 2
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**