

Anti-Uroplakin III antibody (STJ16101371)

STJ16101371

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Uroplakin III is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:50-1:100)
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7380
Gene Symbol	UPK3A
Uniprot ID	UPK3A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

