

Anti-CD9 antibody (STJ16101389)

STJ16101389

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD9 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry and Immunofluorescence research applications.
Applications	WB, IP, FC, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	928
Gene Symbol	CD9
Uniprot ID	CD9_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	