

## Anti-Mycobacterium tuberculosis antibody (STJ16101390)

STJ16101390

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Mycobacterium tuberculosis is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	IHC-P, IHC-F
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Bacteria

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.1% sodium azide
<b>Isotype</b>	
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>