

Anti-c-MET (Hepatocyte Growth Factor Receptor) antibody (STJ16101402)

STJ16101402

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-c-MET (Hepatocyte Growth Factor Receptor) is suitable for use in Immunohistochemistry, Immunocytochemistry and Immunofluorescence research applications.
Applications	IHC-P, ICC, IF
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P ICC IF (1:250)
Formulation	PBS with 1% BSA and 0.1% sodium azide
Isotype	IG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

