

Anti-CD120b antibody (STJ16101405)

STJ16101405

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TNF-RII p75/p80 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IF, IP, WB
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	