

## Anti-ABCC2 antibody (STJ16101427)

STJ16101427

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-MRP2 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
<b>Applications</b>	ICC, IHC-F, IHC-P, FC, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ICC IHC-F (1:20) IHC-P FC WB (1:20-50)
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>