

Anti-Alpha-Synuclein antibody (STJ16101679)

STJ16101679

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Alpha-Synuclein is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC-P |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | IHC-P (1:20-1:40) |
| Formulation | Tissue culture supernatant with 15mM sodium azide |
| Isotype | IgG1k |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 6622 |
| Gene Symbol | SNCA |
| Uniprot ID | SYUA_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |

