

Anti-Beta-Dystroglycan antibody (STJ16101683)

STJ16101683

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta-Dystroglycan is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human, Mouse, Rat, Rabbit, Dog, Chicken, Hamster, Toad

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Lyophilized
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081