

Anti-Caveolin-1 antibody (STJ16101692)

STJ16101692

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Caveolin-1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:200) WB
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	857
Gene Symbol	CAV1
Uniprot ID	CAV1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

