

Anti-CD123 antibody (STJ16101696)

STJ16101696

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD123 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:100)
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	8872
Gene Symbol	CDC123
Uniprot ID	CD123_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

