

Anti-Cyclin D1 antibody (STJ16101725)

STJ16101725

GENERAL INFORMATION

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Cyclin D1 is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC-P |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------------------|
| Clonality | Monoclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | IHC-P (1:50) |
| Formulation | Tissue culture supernatant with Sodium azide |
| Isotype | IgG2a |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |

