

Anti-Folate Receptor Alpha antibody (STJ16101741)

STJ16101741

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Folate Receptor Alpha is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:100)
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2348
Gene Symbol	FOLR1
Uniprot ID	FOLR1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

