

Anti-Prostatic Acid Phosphatase antibody (STJ16101781)

STJ16101781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Prostatic Acid Phosphatase is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P, IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:100) IHC-F
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	55
Gene Symbol	ACP3
Uniprot ID	PPAP_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

