

Anti-Utrophin antibody (STJ16101818)

STJ16101818

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Utrophin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human, Rat, Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Lyophilized
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7402
Gene Symbol	UTRN
Uniprot ID	UTRN_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	