

## Anti-CD45RO antibody (STJ16101830)

STJ16101830

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-CD45RO is suitable for use in Immunohistochemistry and Immunohistochemistry research applications. |
| <b>Applications</b>      | IHC-P, IHC-F  |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal                                   |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | Unconjugated                                 |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      | IHC-P (1:200-1:800)<br>IHC-F                 |
| <b>Formulation</b>         | Tissue culture supernatant with sodium azide |
| <b>Isotype</b>             | IgG2a  |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt.                 |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            |  |
| <b>Gene Symbol</b>        |  |
| <b>Uniprot ID</b>         |  |
| <b>Immunogen</b>          |  |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |