

Anti-ACTA1 antibody [HHF35] (STJ180000)

STJ180000

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HHF35
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Buffer with protein carrier and preservative
Isotype	IgG1/Kappa
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	58
Gene Symbol	ACTA1
Uniprot ID	ACTS_HUMAN
Immunogen	SDS extract of human myocardium
Immunogen Region	
Specificity	Positive control: Skeletal muscle
Immunogen Sequence	



Human uterus stained with anti-MSA antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of smooth muscle.