

Anti-MDM2 antibody (154-167aa) [SMP14] (STJ180166)

STJ180166

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MDM2 (154-167aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMP14
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:50-100
Formulation	Buffer with protein carrier and preservative
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	4193
Gene Symbol	MDM2
Uniprot ID	MDM2_HUMAN
Immunogen	Amino acids 154-167 of MDM2 of human origin
Immunogen Region	154-167aa
Specificity	Positive control: Liposarcoma
Immunogen Sequence	



Human atypical lipomatous tumor stained with anti-MDM2 antibody using peroxidase-conjugate and DAB chromogen. Note the nuclear staining of tumor cells.