

Anti-TYMS antibody [TS106] (STJ180261)

STJ180261

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TYMS is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TS106
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:50
Formulation	Buffer with protein carrier and preservative
Isotype	IgG1/Kappa
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	7298
Gene Symbol	TYMS
Uniprot ID	TYSY_HUMAN
Immunogen	Recombinant human thymidylate synthase
Immunogen Region	
Specificity	Positive control: 5-FU resistant cancer
Immunogen Sequence	



Human colon carcinoma stained with anti-thymidylate synthase (TS) antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of tumor cells.