

Anti-AMACR+p63 antibody [13H4+ZM70] (STJ180269)

STJ180269

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse, Rabbit monoclonal antibody anti-AMACR+p63 is suitable for use in Immunohistochemistry research applications.
Applications	IHC
Host/Source	Mouse, Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	13H4+ZM70
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:50
Formulation	Buffer with protein carrier and preservative
Isotype	IgG, IgG1a/k
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	23600
Gene Symbol	AMACR
Uniprot ID	AMACR_HUMAN
Immunogen	Human AMACR (P504S) and p63 polypeptides
Region	
Specificity	Positive control: Prostate carcinoma
Immunogen Sequence	



Human prostate carcinoma stained anti-AMACR and anti-p63 antibodies using peroxidase-conjugate and DAB chromogen. Note cytoplasmic and luminal circumferential staining of carcinoma glands with AMACR and basal cell nuclear staining with p63.