

Anti-CD11c antibody (637-827aa) [ZM103] (STJ180337)

STJ180337

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11c (637-827aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM103
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Buffer with protein carrier and preservative
Isotype	IgG1/Kappa
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	3687
Gene Symbol	ITGAX
Uniprot ID	ITAX_HUMAN
Immunogen	Recombinant human ITGAX protein fragment (aa 637-827)
Immunogen Region	637-827aa
Specificity	Positive control: Tonsil, lymph node
Immunogen Sequence	



Human tonsil stained with anti-CD11c antibody using peroxidase-conjugate and DAB chromogen. Note strong membranous staining of lymphocytes.