

Anti-Pan-Cytokeratin antibody [OSCAR] (STJ180390)

STJ180390

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Pan-Cytokeratin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OSCAR
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:25-100
Formulation	Buffer with protein carrier and preservative
Isotype	IgG2a/Kappa
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Crude cytokeratin extract prepared from RT-4 and MCF-7 cells
Immunogen Region	
Specificity	Positive control: Squamous cell carcinoma, prostate
Immunogen Sequence	



Human liver stained with anti-pan cytokeratin OSCAR antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of ductal cells (strong) and hepatocytes (weak).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081