

## Anti-CD19 antibody (96-281aa) [ZR212] (STJ180396)

STJ180396

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-CD19 (96-281aa) is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-p
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	ZR212
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	1:100-200
<b>Formulation</b>	Buffer with protein carrier and preservative
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">930</a>
<b>Gene Symbol</b>	<a href="#">CD19</a>
<b>Uniprot ID</b>	<a href="#">CD19_HUMAN</a>
<b>Immunogen</b>	Recombinant fragment of human CD19 protein (around aa 96-281)
<b>Immunogen Region</b>	96-281aa
<b>Specificity</b>	Positive control: Lymph node or tonsil
<b>Immunogen Sequence</b>	



Human tonsil stained with anti-CD19 antibody using peroxidase-conjugate and DAB chromogen. Note the membranous staining of B-cells.