

Anti-Treponema Pallidum antibody (STJ180400)

STJ180400

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Treponema Pallidum is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Rabbit
Reactivity	T.pallidum

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Buffer with protein carrier and preservative
Isotype	Rabbit IgG
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Treponema pallidum
Immunogen Region	
Specificity	Positive control: Treponema pallidum-infected tissue
Immunogen Sequence	



Human lung infected with treponema organisms stained with anti-treponema pallidum using peroxidase-conjugate and DAB chromogen. Note the clusters of slender organisms.