

Anti-SOX-10 antibody (115-269aa) [ZR275] (STJ180423)

STJ180423

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-SOX-10 (115-269aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR275
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Buffer with protein carrier and preservative
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	6663
Gene Symbol	SOX10
Uniprot ID	SOX10_HUMAN
Immunogen	Recombinant human SOX10 protein fragment (around aa 115-269)
Immunogen Region	115-269aa
Specificity	Positive control: Melanoma
Immunogen Sequence	



Human malignant melanoma stained with anti-SOX-10 antibody using peroxidase-conjugate and DAB chromogen. Note the nuclear staining of tumor cells.