

Anti-ERG antibody (450-C-Term) [ZR331] (STJ180452)

STJ180452

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-ERG (450-C-Term) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR331
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	
Formulation	Buffer with protein carrier and preservative
Isotype	IgG
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	2078
Gene Symbol	ERG
Uniprot ID	ERG_HUMAN
Immunogen	Synthetic peptide corresponding to residues within aa 450 to C-terminus in human ERG protein
Immunogen Region	450-C-Term
Specificity	Positive control: Prostate adenocarcinoma
Immunogen Sequence	



Stained with anti-ERG using peroxidase-conjugate and DAB chromogen. Note nuclear staining of vascular endothelial cells.