

Anti-AHNAK antibody (5720-5800) (STJ190105)

STJ190105

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79026
Gene Symbol	AHNAK
Uniprot ID	AHNK_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 5720-5800
Immunogen Region	5720-5800
Specificity	AHNAK polyclonal antibody (Neuroblast Differentiation-Associated Protein Ahnak) binds to endogenous Neuroblast Differentiation-Associated Protein Ahnak at the amino acid region 5720-5800.
Immunogen Sequence	