

Anti-SMC3 antibody (1010-1090) (STJ190157)

STJ190157

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Structural Maintenance Of Chromosomes Protein 3 (1010-1090) is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB, ELISA |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:500-2000 ELISA 1:5000-20000 |
| Formulation | PBS, 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 9126 |
| Gene Symbol | SMC3 |
| Uniprot ID | SMC3_HUMAN |
| Immunogen | Synthesized peptide derived from human protein at aa range 1010-1090 |
| Immunogen Region | 1010-1090 |
| Specificity | SMC3 polyclonal antibody (Structural Maintenance Of Chromosomes Protein 3) binds to endogenous Structural Maintenance Of Chromosomes Protein 3 at the amino acid region 1010-1090. |
| Immunogen Sequence | |