

Anti-BCL3 antibody (340-420) (STJ190158)

STJ190158

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-B-Cell Lymphoma 3 Protein (340-420) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	602
Gene Symbol	BCL3
Uniprot ID	BCL3_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 340-420
Immunogen Region	340-420
Specificity	BCL3 polyclonal antibody (B-Cell Lymphoma 3 Protein) binds to endogenous B-Cell Lymphoma 3 Protein at the amino acid region 340-420.
Immunogen Sequence	