

## Anti-MAPK4 antibody (410-490) (STJ190195)

STJ190195

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Mitogen-Activated Protein Kinase 4 (410-490) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Formulation</b>	PBS, 50% Glycerol and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	5596
<b>Gene Symbol</b>	MAPK4
<b>Uniprot ID</b>	MK04_HUMAN
<b>Immunogen</b>	Synthesized peptide derived from human protein at aa range 410-490
<b>Immunogen Region</b>	410-490
<b>Specificity</b>	MAPK4 polyclonal antibody (Mitogen-Activated Protein Kinase 4) binds to endogenous Mitogen-Activated Protein Kinase 4 at the amino acid region 410-490.
<b>Immunogen Sequence</b>	