

Anti-TRIM35 antibody (320-400) (STJ190262)

STJ190262

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 Ubiquitin-Protein Ligase Trim35 (320-400) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23087
Gene Symbol	TRIM35
Uniprot ID	TRI35_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 320-400
Immunogen Region	320-400
Specificity	TRIM35 polyclonal antibody (E3 Ubiquitin-Protein Ligase Trim35) binds to endogenous E3 Ubiquitin-Protein Ligase Trim35 at the amino acid region 320-400.
Immunogen Sequence	