

Anti-COL9A2 antibody (540-620) (STJ190272)

STJ190272

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Collagen Alpha-2 (I α Chain (540-620) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1298
Gene Symbol	COL9A2
Uniprot ID	CO9A2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 540-620
Immunogen Region	540-620
Specificity	COL9A2 polyclonal antibody (Collagen Alpha-2 (I α Chain) binds to endogenous Collagen Alpha-2 (I α Chain at the amino acid region 540-620.
Immunogen Sequence	