

Anti-CYP4F3 antibody (400-480) (STJ190287)

STJ190287

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 4f3 (400-480) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4051
Gene Symbol	CYP4F3
Uniprot ID	CP4F3_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 400-480
Immunogen Region	400-480
Specificity	CYP4F3 polyclonal antibody (Cytochrome P450 4f3) binds to endogenous Cytochrome P450 4f3 at the amino acid region 400-480.
Immunogen Sequence	