

Anti-EPS15L1 antibody (350-430) (STJ190299)

STJ190299

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Epidermal Growth Factor Receptor Substrate 15-Like 1 (350-430) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58513
Gene Symbol	EPS15L1
Uniprot ID	EP15R_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 350-430
Immunogen Region	350-430
Specificity	EPS15L1 polyclonal antibody (Epidermal Growth Factor Receptor Substrate 15-Like 1) binds to endogenous Epidermal Growth Factor Receptor Substrate 15-Like 1 at the amino acid region 350-430.
Immunogen Sequence	