

Anti-FAT1 antibody (1140-1220) (STJ190301)

STJ190301

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protocadherin Fat 1 (1140-1220) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2195
Gene Symbol	FAT1
Uniprot ID	FAT1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1140-1220
Immunogen Region	1140-1220
Specificity	FAT1 polyclonal antibody (Protocadherin Fat 1) binds to endogenous Protocadherin Fat 1 at the amino acid region 1140-1220.
Immunogen Sequence	