

Anti-KDM5A antibody (720-800) (STJ190343)

STJ190343

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lysine-Specific Demethylase 5a (720-800) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5927
Gene Symbol	KDM5A
Uniprot ID	KDM5A_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 720-800
Immunogen Region	720-800
Specificity	KDM5A polyclonal antibody (Lysine-Specific Demethylase 5a) binds to endogenous Lysine-Specific Demethylase 5a at the amino acid region 720-800.
Immunogen Sequence	