

Anti-MYLK2 antibody (370-450) (STJ190371)

STJ190371

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Myosin Light Chain Kinase 2-Skeletal/Cardiac Muscle (370-450) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85366
Gene Symbol	MYLK2
Uniprot ID	MYLK2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 370-450
Immunogen Region	370-450
Specificity	MYLK2 polyclonal antibody (Myosin Light Chain Kinase 2-Skeletal/Cardiac Muscle) binds to endogenous Myosin Light Chain Kinase 2-Skeletal/Cardiac Muscle at the amino acid region 370-450.
Immunogen Sequence	