

Anti-CYB5R4 antibody (380-460) (STJ190373)

STJ190373

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome B5 Reductase 4 (380-460) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51167
Gene Symbol	CYB5R4
Uniprot ID	NB5R4_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 380-460
Immunogen Region	380-460
Specificity	CYB5R4 polyclonal antibody (Cytochrome B5 Reductase 4) binds to endogenous Cytochrome B5 Reductase 4 at the amino acid region 380-460.
Immunogen Sequence	