

Anti-PLCD1 antibody (90-170) (STJ190380)

STJ190380

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Delta-1 (90-170) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5333
Gene Symbol	PLCD1
Uniprot ID	PLCD1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 90-170
Immunogen Region	90-170
Specificity	PLCD1 polyclonal antibody (1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Delta-1) binds to endogenous 1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Delta-1 at the amino acid region 90-170.
Immunogen Sequence	