

Anti-STRN4 antibody (250-330) (STJ190418)

STJ190418

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Striatin-4 (250-330) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29888
Gene Symbol	STRN4
Uniprot ID	STRN4_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 250-330
Immunogen Region	250-330
Specificity	STRN4 polyclonal antibody (Striatin-4) binds to endogenous Striatin-4 at the amino acid region 250-330.
Immunogen Sequence	