

Anti-SATB2 antibody (650-730) (STJ190433)

STJ190433

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dna-Binding Protein Satb2 (650-730) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23314
Gene Symbol	SATB2
Uniprot ID	SATB2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 650-730
Immunogen Region	650-730
Specificity	SATB2 polyclonal antibody (Dna-Binding Protein Satb2) binds to endogenous Dna-Binding Protein Satb2 at the amino acid region 650-730.
Immunogen Sequence	