

## Anti-ADAM30 antibody (400-480) (STJ190513)

STJ190513

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Disintegrin And Metalloproteinase Domain-Containing Protein 30 (400-480) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution</b>	WB 1:500-2000
<b>Range</b>	ELISA 1:5000-20000
<b>Formulation</b>	PBS, 50% Glycerol and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">11085</a>
<b>Gene Symbol</b>	<a href="#">ADAM30</a>
<b>Uniprot ID</b>	<a href="#">ADA30_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human protein at aa range 400-480
<b>Immunogen Region</b>	400-480
<b>Specificity</b>	ADAM30 polyclonal antibody (Disintegrin And Metalloproteinase Domain-Containing Protein 30) binds to endogenous Disintegrin And Metalloproteinase Domain-Containing Protein 30 at the amino acid region 400-480.
<b>Immunogen Sequence</b>	