

Anti-ADAMTS14 antibody (190-270) (STJ190518)

STJ190518

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-A Disintegrin And Metalloproteinase With Thrombospondin Motifs 14 (190-270) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	140766
Gene Symbol	ADAMTS14
Uniprot ID	ATS14_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 190-270
Immunogen Region	190-270
Specificity	ADAMTS14 polyclonal antibody (A Disintegrin And Metalloproteinase With Thrombospondin Motifs 14) binds to endogenous A Disintegrin And Metalloproteinase With Thrombospondin Motifs 14 at the amino acid region 190-270.
Immunogen Sequence	