

Anti-AGXT2 antibody (340-420) (STJ190527)

STJ190527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Alanine--Glyoxylate Aminotransferase 2-Mitochondrial (340-420) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64902
Gene Symbol	AGXT2
Uniprot ID	AGT2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 340-420
Immunogen Region	340-420
Specificity	AGXT2 polyclonal antibody (Alanine--Glyoxylate Aminotransferase 2-Mitochondrial) binds to endogenous Alanine--Glyoxylate Aminotransferase 2-Mitochondrial at the amino acid region 340-420.
Immunogen Sequence	